

# CardioVision: Predicting Cardiac Arrest Risk

An Ensemble Approach Utilizing Wearable Physiological Data




# Introduction

**Problem: Sudden Cardiac Arrest**



- Sudden Cardiac Arrest (SCA) often occurs without warning symptoms.
- ~350,000 SCA cases occur annually in the United States.
- Survival rate is below 10% in most cases.
- Early risk prediction is critical for improving survival outcomes.
- Current clinical approaches often lack real-time, proactive detection.

**Motivation: Personal Connection and Wearable Technology**



- A personal loss in our team due to cardiac arrest sparked the project.
- Wearable tech like the Apple Watch collects continuous health data (e.g., heart rate, ECG).
- Opportunities for real-time cardiac risk monitoring.
- Millions of users worldwide.

**Goals: Ensemble Model for Real-Time Prediction**



- Build a machine learning model to predict cardiac arrest risk early.
- Classify risk into three levels: Low, Medium, High.
- Combine outputs from multiple specialized models for more accurate prediction.
- Enable real-time risk scoring using physiological data.
- Integrate seamlessly with Apple Watch for proactive alerts and interventions.

**Team and Mentorship**

**Mentors**



Vineel Nagisetty  
Jaemar Miller  
Varun Garde  
Luhhsan Elankumaran  
Nirhar Dey

**LNV**




# Problem: Sudden Cardiac Arrest



- Sudden Cardiac Arrest (SCA) often occurs without warning symptoms.
- **~350,000 SCA cases occur annually** in the United States.
- Survival rate is **below 10%** in most cases.
- Early risk prediction is critical for improving survival outcomes.
- Current clinical approaches often lack real-time, proactive detection.

# Motivation: Personal Connection and Wearable Technology

- A **personal loss** in our team due to cardiac arrest sparked the project.
- **Wearable tech** like the **Apple Watch** collects continuous health data (e.g., heart rate, ECG).
- Opportunities for **real-time cardiac risk monitoring**
- **Millions of users worldwide**



# Goals: Ensemble Model for Real-Time Prediction



- Build a **machine learning model** to predict cardiac arrest risk early.
- Classify risk into **three levels**: Low, Medium, High.
- Combine outputs from **multiple specialized models** for more accurate prediction.
- Enable **real-time risk scoring** using physiological data.
- Integrate seamlessly with **Apple Watch** for proactive **alerts and interventions**.

LNV

# Team and Mentorship

## Mentors



Vineel Nagisetty



Jaemar Millar



Lukhsaan  
Elankumaran



Nirjhar  
Dey



Varun  
Gande

# CardioVision: Predicting Cardiac Arrest Risk

An Ensemble Approach Utilizing Wearable Physiological Data



# Data Sources

### Apple Healthkit API

- Mock JSON Data generated based on healthkit metrics and ranges:
  - Heart Rate (HR)
  - Resting Heart Rate (RHR)
  - High Heart Rate Events (HHR)
  - Heart Rate Variability (HRV)
  - Electrocardiogram (ECG)




### MIT-BIH Arrhythmia Database

- 48 half-hour ECG recordings from 47 patients sampled at 360 Hz
- File Types:
  - Raw Signal (.dat): Raw ECG waveform data
  - Header File (.hea): Metadata file with details like sampling rate, signal names, and signal gains.
  - Annotation File (.atr): Annotations for each beat (normal, arrhythmic, artifact).
  - X-Wave Annotation (.xws): Extended waveform signal data for advanced signal analysis.



### MIMIC-III Waveform Database

- 67,000 ICU waveform records linked to clinical data
- Includes ECG, ABP, respiration, SpO<sub>2</sub> signals
- 48 hours of ECG recordings sampled at 125 Hz
- File Types:
  - Raw Signal (.dat): Raw ECG waveform data
  - Header File (.hea): Metadata file with details like sampling rate, signal names, and signal gains



### Sudden Cardiac Death Holter Database

- ECG recordings from patients at risk of sudden cardiac death (SCD).
- 24-hour Holter ECG data, sampled at 250 Hz
- File Types:
  - Raw Signal (.dat): Raw ECG waveform data
  - Header File (.hea): Metadata file with details like sampling rate, signal names, and signal gains
  - Annotation File (.atr): Annotations for each beat (normal, arrhythmic, artifact)




### St. Petersburg INCART 12-lead Arrhythmia Database

- 75 half-hour 12-lead ECG recordings from 32 patients sampled at 257 Hz with over 175,000 annotated beats
- Patient conditions include atrial fibrillation, SVT, MI, bundle branch blocks, and ischemic events.
- File Types:
  - Raw Signal (.dat): Raw ECG waveform data
  - Header File (.hea): Metadata file with details like sampling rate, signal names, and signal gains
  - Annotation File (.atr): Annotations for each beat (normal, arrhythmic, artifact)




### Out-of-Hospital Cardiac Arrest Database

- 260 ECG recordings from OHCA patients treated by EMS
- Each record includes 9 sec pre-shock + 1 min post-shock ECG segments.
- Each case labeled as successful (ROEA) or unsuccessful (NoROEA) defibrillation by 3 cardiologists.
- File Types:
  - PDF (.pdf): Scanned ECG waveform images (anonymized, printed, scanned).
  - Text (.txt): Digitized ECG waveforms with amplitude-time data.
  - Excel (.xls): Extracted feature values in tabular format (time, frequency, wavelet, and nonlinear metrics).



# Apple Healthkit API

- Mock **JSON** Data generated based on healthkit metrics and ranges:
  - **Heart Rate (HR)** :
  - **Resting Heart Rate (RHR)**
  - **High Heart Rate Events (HHR)**
  - **Heart Rate Variability (HRV)**
  - **Electrocardiogram (ECG)**



# MIT-BIH Arrhythmia Database

- 48 half-hour ECG recordings from 47 patients sampled at 360 Hz
- File Types:
  - Raw Signal (.dat): Raw ECG waveform data.
  - Header File (.hea): Metadata file with details like sampling rate, signal names, and signal gains.
  - Annotation File (.atr): Annotations for each beat (normal, arrhythmic, artifact).
  - X-Wave Annotation (.xws): Extended waveform signal data for advanced signal analysis.




# MIMIC-III Waveform Database

- 67,000 ICU waveform records linked to clinical data
- Includes ECG, ABP, respiration, SpO<sub>2</sub> signals
- 48 hours of ECG recordings sampled at 125 Hz
- File Types:
  - Raw Signal (.dat): Raw ECG waveform data.
  - Header File (.hea): Metadata file with details like sampling rate, signal names, and signal gains.



# Sudden Cardiac Death Holter Database



- ECG recordings from patients at risk of sudden cardiac death (SCD).
- 24-hour Holter ECG data, sampled at 250 Hz.
- File Types:
  - Raw Signal (.dat): Raw ECG waveform data.
  - Header File (.hea): Metadata file with details like sampling rate, signal names, and signal gains.
  - Annotation File (.atr): Annotations for each beat (normal, arrhythmic, artifact).

# St. Petersburg INCART 12-lead Arrhythmia Database

- 75 half-hour **12-lead ECG** recordings from 32 patients sampled at **257 Hz** with over 175,000 annotated beats
- Patient conditions include **atrial fibrillation, SVT, MI, bundle branch blocks, and ischemic events.**




- **File Types:**
  - **Raw Signal (.dat):** Raw ECG waveform data.
  - **Header File (.hea):** Metadata file with details like sampling rate, signal names, and signal gains.
  - **Annotation File (.atr):** Annotations for each beat (normal, arrhythmic, artifact).

# Out-of-Hospital Cardiac Arrest Database

- 260 ECG recordings from OHCA patients treated by EMS.
- Each record includes 9 sec pre-shock + 1 min post-shock ECG segments.
- Each case labeled as successful (ROEA) or unsuccessful (NoROEA) defibrillation by 3 cardiologists.

- File Types:

- PDF (.pdf): Scanned ECG waveform images (anonymized, printed, scanned).
- Text (.txt): Digitized ECG waveforms with amplitude-time data.
- Excel (.xls): Extracted feature values in tabular format (time, frequency, wavelet, and nonlinear metrics).



# CardioVision: Predicting Cardiac Arrest Risk

An Ensemble Approach Utilizing Wearable Physiological Data



# Methods and Pipeline





## CardioVision

### Initial Risk Prediction

- If Prediction = 0 :
  - No risk, continue monitoring
- If Prediction = 1:
  - Prompt user to record ECG:

### Final Risk Prediction

- If Prediction = 0 :
  - False Alarm
- If Prediction = 1:
  - Cardiac Symptoms, Monitor
- If Prediction = 2:
  - Arrest Risk, Contact EMS



HR, RHR HHR, HRV



FastAPI

HR, RHR HHR, HRV



H

# Initial Risk Model

## Initial HealthKit Metrics

### Datasource: Mock Healthkit Data

- Mock HealthKit data generated as JSON files in two categories: "risk" and "no\_risk"
- Metrics:
  - Heart Rate (HR): 60-100 bpm (no risk), 110-140 bpm (risk).
  - Heart Rate Variability (HRV): 55-100 ms (no risk), 10-40 ms (risk).
  - Resting Heart Rate (RHR): 55-85 bpm (no risk), 90-110 bpm (risk).
  - High Heart Rate Events (HHR): 0 (no risk), 1-5 events (risk).
- Ranges are based on established values from the [American Heart Association](#), [Kubios HRV](#), [WebMD](#), and [Harvard Health](#)

### Model: Random Forest Classifier

- **Random Forest:** Aggregates predictions from multiple decision trees for robust classification.
- **Decision Trees:** Individual tree classifiers that split data based on conditions.
- 100 Decision Trees

### Results

Mock Healthkit Data

	Normal	Arrhythmia
Normal	100	0
Arrhythmia	0	100

### Training:

- Loaded **JSON** data into a **DataFrame** with 4 features (HR, HRV, RHR, HHR).
- Split into training (80%) and testing (20%) sets.
- Trained **Random Forest model** on the training set.

# Initial HealthKit Metrics

## Datasource: Mock Healthkit Data

- Mock HealthKit data generated as JSON files in two categories: "risk" and "no\_risk"
- Metrics:
  - **Heart Rate (HR)**: 60-100 bpm (**no risk**), 110-140 bpm (**risk**).
  - **Heart Rate Variability (HRV)**: 55-100 ms (**no risk**), 10-40 ms (**risk**).
  - **Resting Heart Rate (RHR)**: 55-85 bpm (**no risk**), 90-110 bpm (**risk**).
  - **High Heart Rate Events (HHR)**: 0 (**no risk**), 1-5 events (**risk**).
- Ranges are based on established values from the **American Heart Association**, **Kubios HRV**, **WebMD**, and **Harvard Health**

## Model: Random Forest Classifier

- **Random Forest**: Aggregates predictions from multiple decision trees for robust classification.
- **Decision Trees**: Individual tree classifiers that split data based on conditions.
- 100 Decision Trees

## Results

Mock Healthkit Data

	Normal	Arrhythmia
Normal	100	0
Arrhythmia	0	100

## Training:

- Loaded **JSON** data into a **DataFrame** with 4 features (**HR**, **HRV**, **RHR**, **HHR**).
- Split into training (80%) and testing (20%) sets.
- Trained **Random Forest model** on the training set.



*Prediction*



# FastAPI

Prediction



# Initial Risk Prediction

- If Prediction = 0 :
  - No risk, continue monitoring
- If Prediction = 1:
  - Prompt user to record ECG:

W ECG







# FastAPI



**ECG**

# Final Risk Model




# Heart Rate (HR)

## Datasource: MIT-BIH

- HR Computation: Calculated as  $HR = 60/RR$  : where RR is the interval between consecutive R-peaks (seconds).

## Model: Rule-Based Threshold

- Binary Classification: **High HR (1) if  $HR > 110$  bpm, Normal HR (0) otherwise.**
- No training required as it is a fixed threshold



## Results

### MIT-BIH

	Normal	Arrhythmia
Normal	23063	0
Arrhythmia	0	12

### INCART

	Normal	Arrhythmia
Normal	58951	0
Arrhythmia	0	78

# High Heart Rate Events (HHR)

## Datasource: MITBIH

- RR intervals calculated from R-peaks in ECG, converted to HR (BPM).

- Features:

- Duration above 150 BPM
- sustained high HR (all > 150 BPM),
- max, min, avg HR,
- HR slope
- spike frequency
- HR standard deviation.

## Model: Random Forest Classifier

- Random Forest: Aggregates predictions from multiple decision trees for robust classification.
- Decision Trees: Individual tree classifiers that split data based on conditions.
- 300 Trees, max depth of 10, balanced class weights

## Results

MIT-BIH

	Normal	Arrhythmia
Normal	59583	0
Arrhythmia	0	53117

INCART

	Normal	Arrhythmia
Normal	20674	0
Arrhythmia	0	14656

## Training:

- Extracted HR features scaled using StandardScaler, balanced with SMOTE (Synthetic Minority Over-sampling - interpolating between minority class samples to introduce variation) and mapped to labels:
  - High HR event = High duration above threshold + sustained High HR
  - Normal event = Lower values or mixed patterns

# Resting Heart Rate (RHR)

## Results

### MIMIC-III

	Normal	Arrhythmia
Normal	96	0
Arrhythmia	0	218

### INCART

	Normal	Arrhythmia
Normal	957	1
Arrhythmia	0	1187

### Datasource: MIMIC-III

- ECG signals processed to detect RR intervals using peak detection.
- Mean HR and standard deviation (filtered for low variation, std < 5 for MIMIC, std < 15 for INCART).


### Model: Logistic Regression

- Linear model that predicts binary outcomes using the sigmoid function

### Training:

- Extracted RHR values standardized with StandardScaler, labeled as high (1) or normal (0) based on a 75 BPM threshold, and split into training/test sets.

$$\frac{e^{(\beta_0 + \beta_1 x)}}{1 + e^{(\beta_0 + \beta_1 x)}}$$



# Heart Rate Variability (HRV)

## Datasource: MIT-BIH

- RR intervals computed from ECG signals
- Calculated 14 advanced HRV features for each 30-beat window:
  - Time-Domain:
    - RMSSD (Root Mean Square of Successive Differences)
    - SDNN (Standard Deviation of NN (Normal-Normal) Intervals)
    - NN50 (Number of consecutive RR intervals differing by more than 50ms,
    - pNN50 (Percentage of NN50 relative to total intervals)
    - Mean, Min, Max, Skew abd Kurtosis of RR Intervals
  - Frequency-Domain:
    - LF (Low-Frequency), HF (High-Frequency), LF/HF ratio using Welch's method.
  - Nonlinear: SD1, SD2 (Short and Long term HRV), Shannon entropy (unpredictability)

## Training:

- HRV features scaled with StandardScaler.
- Models trained using optimized hyperparameters (RandomizedSearchCV).
- Ensemble model trained with StackingClassifier using predictions from XGBoost, LightGBM, and CatBoost.

## Results

### MIT-BIH

	Normal	Arrhythmia
Normal	18735	12491
Arrhythmia	2050	2298

### INCART

	Normal	Arrhythmia
Normal	25135	37154
Arrhythmia	4156	8895

Precision	0.87	0.22
Recall	0.60	0.53

Precision	0.90	0.15
Recall	0.15	0.68

Accuracy	0.62
----------	------

Accuracy	0.46
----------	------

## Models:

- Gradient Boosting: Combines weak models (trees) to improve accuracy.
  - XGBoost: Fast, depth-controlled trees.
  - LightGBM: Leaf-wise growth, fast training.
  - CatBoost: Categorical feature support.
- Ensemble (Stacking Classifier):
  - Combines predictions from XGBoost, LightGBM, and CatBoost using Logistic Regression for improved accuracy.

# Datasource: MITBIH/Holter

- ECG beat segments from raw signals classified into:
  - Low Risk: Normal rhythms:
  - N (Normal)
  - L (Left bundle branch block)
  - R (Right bundle branch block)
  - e (Atrial escape)
  - j (Nodal escape).
- Medium Risk: Supraventricular arrhythmias
  - A (Atrial premature)
  - S (Supraventricular)
  - a (Aberrated atrial)
  - J (Nodal premature)
  - ? (Unclassifiable).
- High Risk: Ventricular arrhythmias:
  - V (Premature ventricular)
  - F (Fusion of ventricular and normal)
  - E (Ventricular escape).

# Electrocardiogram (ECG)

## Results

### OHCA

	Low	Medium	High
Low	0	0	0
Medium	0	0	0
High	80	39	3455

Precision	0.00	0.00	0.99
Recall	0.00	0.00	0.97
Accuracy			0.97

### MIT-BIH + Holter

	Low	Medium	High
Low	24177	75	254
Medium	0	475	1
High	170	104	23650

Precision	0.99	0.73	0.99
Recall	0.99	0.99	0.99
Accuracy			0.99

### INCART

	Low	Medium	High
Low	91664	43100	18832
Medium	545	1355	59
High	3473	1001	15751

Precision	0.96	0.03	0.45
Recall	0.60	0.69	0.78
Accuracy			0.62

## Model: BiLSTM

- Bidirectional Long Short-Term Memory (BiLSTM) - a type of recurrent neural network (RNN) that learns from both past and future contexts in sequential data, enhancing pattern recognition - with 3 hidden layers (128 hidden units each), followed by dropout (0.3) and a fully connected output layer.
- Captures bidirectional temporal dependencies in ECG waveforms for robust pattern recognition.
- Takes 250-sample ECG beat segments as input and outputs multiclass classification (0: Low Risk, 1: Medium Risk, 2: High Risk).

## Training:

- CrossEntropyLoss with class weights ensures all risk levels are learned equally
- Adam optimizer with a learning rate of 0.0005 for efficient training.
- Medium-risk beats are augmented 10x with random transformations to ensure balanced training.
- 30 epochs (full training passes through data), monitoring validation accuracy to prevent overfitting
- Trained in small groups of 32 beats for faster, stable learning

# Final Risk Model

## Heart Rate (HR)

- Datasource:** MIT-BIH  
 - HR Computation: Calculated as  $HR = 60/RR$  where  $RR$  is the interval between consecutive R-peaks (seconds).  
 - Binary Classification: High HR (1) if  $HR > 110$  bpm, Normal HR (0) otherwise.  
 - No training required as it is a fixed threshold.

**Model: Rule-Based Threshold**

# Final Ensemble

## Extraction Process:

- Extracts **HR**, **HRV**, **RHR**, **HHR**, and **ECG** features using **BiLSTM** (3 layers, 128 hidden units, bidirectional).
- Applies **SMOTE** (Synthetic Minority Oversampling) for balanced training data.

## Model: Adaptive Meta-Learner

- Combines predictions from **HR**, **HRV**, **RHR**, **HHR**, and **ECG** (**BiLSTM**) using a **Meta-Learner** (**Neural Network**).
- False Negative Reduction:** Adaptive voting adjusts high-risk predictions using Focal Loss.

## Results

OHCA

	Low	Medium	High
Low	0	0	0
Medium	0	0	0
High	31	19	3524

MIT-BIH + Holter

	Low	Medium	High
Low	24200	60	246
Medium	0	475	1
High	62	38	23824

INCART

	Low	Medium	High
Low	105679	29417	18500
Medium	500	1400	59
High	2499	1001	16725

	Precision	Recall	Accuracy
Precision	0.00	0.00	0.99
Recall	0.00	0.00	0.98

	Precision	Recall	Accuracy
Precision	0.99	0.89	0.99
Recall	0.99	0.99	0.99

	Precision	Recall	Accuracy
Precision	0.96	0.05	0.47
Recall	0.69	0.71	0.83

## Training:

- Fine-tunes **BiLSTM** with feedback samples (high-risk FN corrections) for **25 epochs** (LR: **0.0001**).
- Optimized with **Focal Loss** ( $\alpha=0.95$ ,  $\gamma=2.0$ ) to further minimize **false negatives**.
- Tracks accuracy, confusion matrix, and false negative rates.



Prediction



# FastAPI




Prediction

# Final Risk Prediction

- If Prediction = 0 :
  - False Alarm
- If Prediction = 1:
  - Cardiac Symptoms, Monitor
- If Prediction = 2:
  - Arrest Risk, Contact EMS

# Demo



● ● ● Series 9 (45mm) – watchOS 10.5



11:31

## CardioVision

Latest Heart Rate

**0 bpm**

Prediction Result

**No Risk**



# CardioVision: Predicting Cardiac Arrest Risk

An Ensemble Approach Utilizing Wearable Physiological Data



## Discussion and Future Directions

### Enhanced Initial Risk Model

- Expanded Metrics: Integrate additional Apple Watch metrics to refine initial risk assessment such as:
  - Heart Rate
  - VO2 Max
  - Blood Oxygen
  - SpO2
  - Steps
- Trend Risk and Resilience for personal insights over time to identify trends and make more informed predictions.



### Improved HRV Model and Dataset

- HRV model struggles with distinguishing arrhythmias due to limited arrhythmia specific data
- Use the PhysioNet MIT-BIH Arrhythmia Database, which is specifically designed for arrhythmia detection



### References

- Apple Watch Data Integration. Apple Developer. Available at: [https://developer.apple.com/documentation/watchkit/integrating\\_your\\_app\\_with\\_the\\_watchkit\\_framework/integrating\\_your\\_app\\_with\\_the\\_watchkit\\_framework](https://developer.apple.com/documentation/watchkit/integrating_your_app_with_the_watchkit_framework/integrating_your_app_with_the_watchkit_framework)
- mHealth Toolkit. This is a Research Toolkit of Evidence-Based Tools for mHealth Applications. Available at: <https://www.mhealthtoolkit.org/>
- Arrhythmia Detection. Wikipedia. Available at: <https://en.wikipedia.org/w/index.php?title=Arrhythmia&oldid=98000000>
- PhysioNet. PhysioNet: Computing in Cardiology. Available at: <https://physionet.org/>
- MIT-BIH Arrhythmia Database. MIT-BIH Arrhythmia Database Version 1.0. PhysioNet. Available at: <https://physionet.org/resource/physionet/1.0.0/arrhythmia-database-12lead>
- American Heart Association. What is Heart Rate? Retrieved last update 2023-05-16 from <https://www.heart.org/en/heart-health/basics/what-is-heart-rate#:~:text=Heart%20rate%20is%20the%20number%20of%20times%20a%20person's%20heart%20beats%20in%20a%20one-minute%20period.>
- American Heart Association. What is Heart Rate Range? Retrieved last update 2023-05-16 from <https://www.heart.org/en/heart-health/basics/what-is-heart-rate#:~:text=Normal%20adult%20heart%20rate%20is%20between%2060%20and%20100%20beats%20per%20minute.>
- American Heart Association. What is Heart Rate Recovery? Retrieved last update 2023-05-16 from <https://www.heart.org/en/heart-health/basics/what-is-heart-rate#:~:text=Heart%20rate%20recovery%20is%20the%20time%20it%20takes%20for%20your%20heart%20rate%20to%20return%20to%20a%20resting%20level%20after%20you%20stop%20exercising.>
- American Heart Association. What is Heart Rate Variability? Retrieved last update 2023-05-16 from <https://www.heart.org/en/heart-health/basics/what-is-heart-rate#:~:text=Heart%20rate%20variability%20is%20the%20amount%20of%20change%20in%20your%20heart%20rate%20over%20a%20short%20period%20of%20time.>

### Advanced Watch App Integration

- Deployment on Apple Watch
- Real-time remote data storage
- Real-time data streaming with Kafka and FastAPI server
- Have user get notification if initial risk identified




### Conclusion: Promising Ensemble Approach

- Accurate Risk Prediction
- Real-Time Assessment without ECG
- Modular Design: Each model independently trainable and replaceable
- Future Work: Further model training and optimization, Apple Watch and other smart watch integration



# Enhanced Initial Risk Model

- **Expanded Metrics:** Integrate additional Apple Watch metrics to refine initial risk assessment such as:
  - HRV
  - VO2 Max
  - Blood Oxygen
  - Respiratory Rate
  - Stress
- Design to track and account for previous inputs over time to identify trends and make more informed predictions



# Advanced Watch App Integration


---


- Deployment on **Apple Watch**
- Reliable and private data storage
- Real-time data streaming with **Kafka** and **FastAPI** server
- Have user get notification if **initial risk** identified



# Improved HRV Model and Dataset

- **HRV model** struggles with distinguishing arrhythmias due to **limited arrhythmia-specific data**.
- Use the **PhysioNet MIT-BIH Atrial Fibrillation Database**, which is specifically designed for arrhythmia detection.





# Conclusion: Promising Ensemble Approach

- Accurate Risk Prediction
- **Real-Time** Assessment without **ECG**
- **Modular Design:** Each model **independently trainable** and **replaceable**
- **Future Directions:** Further model training and optimization, Apple Watch and other smart watch integration

# References

- Apple HealthKit Documentation: Apple Developer. Available at: <https://developer.apple.com/documentation/healthkit>
- OHCA Database: “Out-of-Hospital Cardiac Arrest Registry,” PMC. Available at: <https://PMC7753135/>
- INCART Database: “The St. Petersburg Institute of Cardiological Technics 12-lead Arrhythmia Database,” PhysioNet. Available at: <https://archive.physionet.org/physiobank/database/incartdb/>
- Holter Database (SDDB): “Sudden Cardiac Death Holter Database,” PhysioNet. Available at: <https://archive.physionet.org/physiobank/database/sddb/>
- MIMIC-III Waveform Database: “MIMIC-III Waveform Database (Version 1.0),” PhysioNet. Available at: <https://physionet.org/content/mimic3fdb/1.0/>
- MIT-BIH Arrhythmia Database: “MIT-BIH Arrhythmia Database (Version 1.0.0),” PhysioNet. Available at: <https://www.physionet.org/content/mitdb/1.0.0/>
- American Heart Association. (n.d.). Tachycardia (Fast Heart Rate). Retrieved from <https://www.heart.org/en/health-topics/arrhythmia/about-arrhythmia/tachycardia--fast-heart-rate>
- Kubios. (n.d.). Heart Rate Variability (HRV) Normal Range. Retrieved from <https://www.kubios.com/blog/heart-rate-variability-normal-range/>
- WebMD. (n.d.). What Is Heart Rate Variability?. Retrieved from <https://www.webmd.com/heart/what-is-heart-rate-variability>
- American Heart Association. (n.d.). Tachycardia (Fast Heart Rate). Retrieved from <https://www.heart.org/en/health-topics/arrhythmia/about-arrhythmia/tachycardia--fast-heart-rate>
- American Heart Association. (n.d.). Tachycardia (Fast Heart Rate). Retrieved from <https://www.heart.org/en/health-topics/arrhythmia/about-arrhythmia/tachycardia--fast-heart-rate>
- Harvard Health Publishing. (2019, March). What your heart rate is telling you. Retrieved from <https://www.health.harvard.edu/heart-health/what-your-heart-rate-is-telling-you>

# Questions & Answers

# CardioVision: Predicting Cardiac Arrest Risk

An Ensemble Approach Utilizing Wearable Physiological Data

